

Product Datasheet

KIF5A Antibody (orb1243466)

Description

KIF5A Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human KIF5A.

Target

KIF5A

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store KIF5A antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

KIF5A antibody can be used for detection of KIF5A by ELISA at 1:12500. KIF5A antibody can be used for detection of KIF5A by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

117 kDa

Uniprot ID

[Q12840](#)

NCBI


[NP_004975](#)

Dilution Range

KIF5A antibody can be used for detection of KIF5A by ELISA at 1:12500. KIF5A antibody can be used for detection of KIF5A by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in IHC on Human Lung.



Antibody used in WB on Human brain at 0....