

Product Datasheet

ACAT2 Antibody (orb1243480)

Description

ACAT2 Antibody


Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human A1BG.
Target	ACAT2
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store A1BG antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	A1BG antibody can be used for detection of A1BG by ELISA at 1:1562500. A1BG antibody can be used for detection of A1BG by western blot at 0.1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	54 kDa
Uniprot ID	P04217
NCBI	NP_570602
Dilution Range	A1BG antibody can be used for detection of A1BG by ELISA at 1:1562500. A1BG antibody can be used for detection of A1BG by western blot at 0.1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in WB on Human Fetal Liver...



Antibody used in WB on Human HepG2 at 1....



Antibody used in WB on Human at 5 ug/ml.