

Product Datasheet

KIF3B Antibody (orb1243489)

Description

KIF3B Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human PARP11.

Target

KIF3B

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store PARP11 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

PARP11 antibody can be used for detection of PARP11 by ELISA at 1:312500. PARP11 antibody can be used for detection of PARP11 by western blot at 1.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

39 kDa

Uniprot ID

[Q9NR21](#)

NCBI

[NP_065100](#)

Dilution Range

PARP11 antibody can be used for detection of PARP11 by ELISA at 1:312500. PARP11 antibody can be used for detection of PARP11 by western blot at 1.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

