

Product Datasheet

CLDN13 Antibody (orb1243527)

Description

CLDN13 Antibody

Species/Host

Rabbit

Reactivity

Mouse

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Human
Jurkat 1 ug...

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of mouse CLDN13.

Target

CLDN13

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store CLDN13 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

CLDN13 antibody can be used for detection of CLDN13 by ELISA at 1:1562500. CLDN13 antibody can be used for detection of CLDN13 by western blot at 1.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

23 kDa

Uniprot ID

[Q9Z0S4](#)

NCBI

[NP_065250](#)

Dilution Range

CLDN13 antibody can be used for detection of CLDN13 by ELISA at 1:1562500. CLDN13 antibody can be used for detection of CLDN13 by western blot at 1.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.