

## Product Datasheet

### ZFP36L1 Antibody (orb1243584)

## Description

ZFP36L1 Antibody

### Species/Host

Rabbit

### Reactivity

Canine, Human

### Conjugation

Unconjugated

### Tested

ELISA, WB

### Applications

Antibody  
used in WB  
on Human  
Placenta at...

### Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ZFP36L1.

### Target

ZFP36L1

### Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

### Form/Appearance

Liquid

### Concentration

batch dependent

### Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store ZFP36L1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

### Note

For research use only

### Application notes

ZFP36L1 antibody can be used for detection of ZFP36L1 by ELISA at 1:62500. ZFP36L1 antibody can be used for detection of ZFP36L1 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

### Clonality

Polyclonal

### MW

36 kDa

### Uniprot ID

[Q07352](#)

### NCBI

[NP\\_004917](#)

### Dilution Range

ZFP36L1 antibody can be used for detection of ZFP36L1 by ELISA at 1:62500. ZFP36L1 antibody can be used for detection of ZFP36L1 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

### Expiration Date

12 months from date of receipt.