

Product Datasheet

FEZF2 Antibody (orb1243634)

Description

FEZF2 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody produced in rabbits immunized with a synthetic peptide corresponding to a region of human ZNF312.

Antibody used in WB on Human Jurkat 2.5 ...

Target

FEZF2

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store ZNF312 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

ZNF312 antibody can be used for detection of ZNF312 by ELISA at 1:312500. ZNF312 antibody can be used for detection of ZNF312 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

50 kDa

Uniprot ID

[Q8TBJ5](#)

NCBI

[NP_060478](#)

Dilution Range

ZNF312 antibody can be used for detection of ZNF312 by ELISA at 1:312500. ZNF312 antibody can be used for detection of ZNF312 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

