

Product Datasheet

UBE2N Antibody (orb1243635)

Description

UBE2N Antibody

Species/Host

Rabbit

Reactivity

Canine, Drosophila, Human, Mouse, Rat, Zebrafish

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human UBE2N.

Target

UBE2N

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store UBE2N antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

UBE2N antibody can be used for detection of UBE2N by ELISA at 1:1562500. UBE2N antibody can be used for detection of UBE2N by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

17 kDa

Uniprot ID

[P61088](#)

NCBI

[NP_003339](#)

Dilution Range

UBE2N antibody can be used for detection of UBE2N by ELISA at 1:1562500. UBE2N antibody can be used for detection of UBE2N by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date


12 months from date of receipt.

110kDa
80kDa

38kDa
21kDa

Antibody used in WB on Human Jurkat cell...

UBE2N



Antibody-RED
DAPI-Blue

See H2-1 Data and Customer Feedback for more information

Antibody used in IHC on Mouse brainstem ...

UBE2N

1 2
250kDa
200kDa
150kDa
100kDa
75kDa
50kDa
35kDa
25kDa
15kDa
10kDa
7.5kDa
5kDa

See H2-1 Data and Customer Feedback for more information

Antibody used in WB on Human NT-2, mouse...