

Product Datasheet

IVNS1ABP Antibody (orb1243640)

Description

IVNS1ABP Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human IVNS1ABP.

Target

IVNS1ABP

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store IVNS1ABP antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

IVNS1ABP antibody can be used for detection of IVNS1ABP by ELISA at 1:312500. IVNS1ABP antibody can be used for detection of IVNS1ABP by western blot at 2.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

72 kDa

Uniprot ID

[Q9Y6Y0](#)

NCBI


[NP_006460](#)

Dilution Range

IVNS1ABP antibody can be used for detection of IVNS1ABP by ELISA at 1:312500. IVNS1ABP antibody can be used for detection of IVNS1ABP by western blot at 2.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in IHC on Human Stomach.



Antibody used in WB on Human HepG2 at 2....