


Product Datasheet

SOX18 Antibody (orb1243651)

Description

SOX18 Antibody


Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SOX18.
Target	SOX18
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store SOX18 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	SOX18 antibody can be used for detection of SOX18 by ELISA at 1:62500. SOX18 antibody can be used for detection of SOX18 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	41 kDa
Uniprot ID	P35713
NCBI	NP_060889
Dilution Range	SOX18 antibody can be used for detection of SOX18 by ELISA at 1:62500. SOX18 antibody can be used for detection of SOX18 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in IHC on Human kidney at ...



Antibody used in WB on Human heart at 0....



Antibody used in IF on Human lymphatic e...