

Product Datasheet

PPP1R13L Antibody (orb1243708)


Description

PPP1R13L Antibody

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human PPP1R13L.
Target	PPP1R13L
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store PPP1R13L antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	PPP1R13L antibody can be used for detection of PPP1R13L by ELISA at 1:62500. PPP1R13L antibody can be used for detection of PPP1R13L by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	89 kDa
Uniprot ID	Q8WUF5
NCBI	NP_006654
Dilution Range	PPP1R13L antibody can be used for detection of PPP1R13L by ELISA at 1:62500. PPP1R13L antibody can be used for detection of PPP1R13L by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.




Antibody used in WB on Transfected 293T ...



See Immunoblot Data and Customer Feedback for more information

Antibody used in WB on Human cell lines ...



See IHC Data and Customer Feedback for more information

Antibody used in IHC on Human A549 at 1:...