



Product Datasheet

PDK4 Antibody (orb1243719)

Description	PDK4 Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC, WB	Antibody used in IHC on Human Muscle.
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human PDK4.	
Target	PDK4	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	Antibody used in WB on Human 293T at 0.2...
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store PDK4 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	PDK4 antibody can be used for detection of PDK4 by ELISA at 1:312500. PDK4 antibody can be used for detection of PDK4 by western blot at 1.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	46 kDa	
Uniprot ID	Q16654	Antibody used in WB on Human Huh7, HepG2...
NCBI	NP_002603	
Dilution Range	PDK4 antibody can be used for detection of PDK4 by ELISA at 1:312500. PDK4 antibody can be used for detection of PDK4 by western blot at 1.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	