

## Product Datasheet

### OLIG2 Antibody (orb1243735)

## Description

OLIG2 Antibody

### Species/Host

Rabbit

### Reactivity

Human, Mouse, Rat

### Conjugation

Unconjugated

### Tested

ELISA, WB

### Applications

Antibody  
used in WB  
on Human  
HepG2 at 0....

### Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human OLIG2.

### Target

OLIG2

### Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

### Form/Appearance

Liquid

### Concentration

batch dependent

### Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store OLIG2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Antibody  
used in IHC  
on Human  
optic nerv...

### Note

For research use only

### Application notes

OLIG2 antibody can be used for detection of OLIG2 by ELISA at 1:312500. OLIG2 antibody can be used for detection of OLIG2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

### Clonality

Polyclonal

### MW

32 kDa

### Uniprot ID

[Q13516](#)

### NCBI

[NP\\_005797](#)

### Dilution Range

OLIG2 antibody can be used for detection of OLIG2 by ELISA at 1:312500. OLIG2 antibody can be used for detection of OLIG2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

### Expiration Date

12 months from date of receipt.