


Product Datasheet

NFIX Antibody (orb1243747)


Description

NFIX Antibody

Species/Host	Rabbit
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human NFIX.
Target	NFIX
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store NFIX antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	NFIX antibody can be used for detection of NFIX by ELISA at 1:312500. NFIX antibody can be used for detection of NFIX by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	49 kDa
Uniprot ID	Q14938
NCBI	NP_002492
Dilution Range	NFIX antibody can be used for detection of NFIX by ELISA at 1:312500. NFIX antibody can be used for detection of NFIX by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in IHC on mouse hippocampus...



Antibody used in IHC on Human Prostate a...



Antibody used in WB on Human Jurkat 0.2-...