

Product Datasheet

FOXA1 Antibody (orb1243773)

Description

FOXA1 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human FOXA1.

Target

FOXA1

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store FOXA1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

FOXA1 antibody can be used for detection of FOXA1 by ELISA at 1:312500. FOXA1 antibody can be used for detection of FOXA1 by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

49 kDa

Uniprot ID

[Q0JTB6](#)

NCBI

[NP_004487](#)

Dilution Range


FOXA1 antibody can be used for detection of FOXA1 by ELISA at 1:312500. FOXA1 antibody can be used for detection of FOXA1 by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in WB on Transfected 293T ...



Antibody used in WB on MCF7 at 1:1000 (L...