

Product Datasheet

HMX1 Antibody (orb1243774)

Description

HMX1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Human
Daudi at 5....

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human HMX1.

Target

HMX1

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store HMX1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Antibody
used in WB
on Human
Daudi at 4....

Note

For research use only

Application notes

HMX1 antibody can be used for detection of HMX1 by ELISA at 1:1562500. HMX1 antibody can be used for detection of HMX1 by western blot at 5.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

39 kDa

Uniprot ID

[Q9NP08](#)

NCBI

[NP_061815](#)

Dilution Range

HMX1 antibody can be used for detection of HMX1 by ELISA at 1:1562500. HMX1 antibody can be used for detection of HMX1 by western blot at 5.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

90 kDa

42 kDa

32 kDa

23 kDa

90 kDa

60 kDa

42 kDa

32 kDa

23 kDa