

Product Datasheet

HMX1 Antibody (orb1243775)

Description

HMX1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human HMX1.

Target

HMX1

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store HMX1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

HMX1 antibody can be used for detection of HMX1 by ELISA at 1:1562500. HMX1 antibody can be used for detection of HMX1 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

39 kDa

Uniprot ID

[Q9NP08](#)

NCBI

[NP_061815](#)

Dilution Range

HMX1 antibody can be used for detection of HMX1 by ELISA at 1:1562500. HMX1 antibody can be used for detection of HMX1 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

90 kDa

42 kDa

32 kDa

23 kDa

Antibody used in WB on Human Daudi at 5....

90 kDa

60 kDa

42 kDa

32 kDa

23 kDa

Antibody used in WB on Human Daudi at 4....



Antibody used in WB on Human cell lines ...