

Product Datasheet

HDAC2 Antibody (orb1243777)

Description

HDAC2 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human HDAC2.

Target

HDAC2

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store HDAC2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

HDAC2 antibody can be used for detection of HDAC2 by ELISA at 1:12500. HDAC2 antibody can be used for detection of HDAC2 by western blot at 0.2-2.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

55 kDa

Uniprot ID

[Q92769](#)

NCBI


[NP_001518](#)

Dilution Range


HDAC2 antibody can be used for detection of HDAC2 by ELISA at 1:12500. HDAC2 antibody can be used for detection of HDAC2 by western blot at 0.2-2.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in IHC on Human Spermatoph...



Antibody used in WB on Human Jurkat 0.2-...