

Product Datasheet

ABT1 Antibody (orb1243781)

Description

ABT1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ABT1.

Target

ABT1

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store ABT1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

ABT1 antibody can be used for detection of ABT1 by ELISA at 1:1562500. ABT1 antibody can be used for detection of ABT1 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

31 kDa

Uniprot ID

[Q9ULW3](#)

NCBI


[NP_037507](#)

Dilution Range


ABT1 antibody can be used for detection of ABT1 by ELISA at 1:1562500. ABT1 antibody can be used for detection of ABT1 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in IHC on Human kidney.



Antibody used in WB on Human heart at 1....