

Product Datasheet

GRM6 Antibody (orb1243782)


Description

GRM6 Antibody

Species/Host	Rabbit
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human GRM6.
Target	GRM6
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store GRM6 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	GRM6 antibody can be used for detection of GRM6 by ELISA at 1:62500. GRM6 antibody can be used for detection of GRM6 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	95 kDa
Uniprot ID	O15303
NCBI	NP_000834
Dilution Range	GRM6 antibody can be used for detection of GRM6 by ELISA at 1:62500. GRM6 antibody can be used for detection of GRM6 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.




Antibody used in WB on Human Brain at 0....



See Immunoblot Data and Customer Feedback for more information

Antibody used in WB on Mouse retina at 1...



Antibody used in WB on Human HepG2 at 1 ...