

Product Datasheet

TAF5L Antibody (orb1243788)

Description	TAF5L Antibody	
Species/Host	Rabbit	
Reactivity	Canine, Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	Antibody used in WB on Human Jurkat 2.5 ...
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human TAF5L.	
Target	TAF5L	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store TAF5L antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	TAF5L antibody can be used for detection of TAF5L by ELISA at 1:62500. TAF5L antibody can be used for detection of TAF5L by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	66 kDa, 37 kDa	
Uniprot ID	O75529	
NCBI	NP_055224	
Dilution Range	TAF5L antibody can be used for detection of TAF5L by ELISA at 1:62500. TAF5L antibody can be used for detection of TAF5L by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	