

Product Datasheet

POU2F3 Antibody (orb1243798)

Description

POU2F3 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human POU2F3.

Target

POU2F3

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store POU2F3 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

POU2F3 antibody can be used for detection of POU2F3 by ELISA at 1:1562500. POU2F3 antibody can be used for detection of POU2F3 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

47 kDa

Uniprot ID

[Q9UKI9](#)

NCBI


[NP_055167](#)

Dilution Range


POU2F3 antibody can be used for detection of POU2F3 by ELISA at 1:1562500. POU2F3 antibody can be used for detection of POU2F3 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.




Antibody used in IHC on Human Muscle.



POU2F3
See IHC Data and Customer Feedback for more information

Antibody used in IHC on Mouse tongue at ...



Antibody used in WB on Human HepG2 at 0....