

Product Datasheet

CIC Antibody (orb1243807)

Description

CIC Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human Clon Channel.

Target

CIC

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store Clon Channel antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

Clon Channel antibody can be used for detection of Clon Channel by ELISA at 1:62500. Clon Channel antibody can be used for detection of Clon Channel by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

164 kDa

Uniprot ID

[Q96RK0](#)

NCBI


[NP_055940](#)

Dilution Range

Clon Channel antibody can be used for detection of Clon Channel by ELISA at 1:62500. Clon Channel antibody can be used for detection of Clon Channel by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in IHC on Human Lung.



Antibody used in WB on Human Jurkat 1.25...