


Product Datasheet

SIRT5 Antibody (orb1243839)


Description

SIRT5 Antibody


Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SIRT5.
Target	SIRT5
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store SIRT5 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	SIRT5 antibody can be used for detection of SIRT5 by ELISA at 1:12500. SIRT5 antibody can be used for detection of SIRT5 by western blot at 1.0-2.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	33 kDa
Uniprot ID	Q9NXA8
NCBI	NP_112534
Dilution Range	SIRT5 antibody can be used for detection of SIRT5 by ELISA at 1:12500. SIRT5 antibody can be used for detection of SIRT5 by western blot at 1.0-2.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in IHC on Human Lung.



Antibody used in IHC on Human kidney.



Antibody used in WB on Human Jurkat 1.0-...