


Product Datasheet

MMP19 Antibody (orb1243870)


Description

MMP19 Antibody

Species/Host	Rabbit
Reactivity	Canine, Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human MMP19.
Target	MMP19
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store MMP19 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	MMP19 antibody can be used for detection of MMP19 by ELISA at 1:62500. MMP19 antibody can be used for detection of MMP19 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	56 kDa
Uniprot ID	Q99542
NCBI	NP_002420
Dilution Range	MMP19 antibody can be used for detection of MMP19 by ELISA at 1:62500. MMP19 antibody can be used for detection of MMP19 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in IHC on Human urinary bl...



Antibody used in IHC on Human urinary bl...



Antibody used in WB on Human HepG2 at 0....