

## Product Datasheet

### ETV5 Antibody (orb1243884)

## Description

ETV5 Antibody

### Species/Host

Rabbit

### Reactivity

Canine, Human, Mouse, Rat

### Conjugation

Unconjugated

### Tested

ELISA, WB

### Applications

### Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ETV5.

### Target

ETV5

### Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

### Form/Appearance

Liquid

### Concentration

batch dependent

### Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store ETV5 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

### Note

For research use only

### Application notes

ETV5 antibody can be used for detection of ETV5 by ELISA at 1:62500. ETV5 antibody can be used for detection of ETV5 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

### Clonality

Polyclonal

### MW

58 kDa, 58 kDa

### Uniprot ID

[P41161](#)

### NCBI


[NP\\_004445](#)

### Dilution Range


ETV5 antibody can be used for detection of ETV5 by ELISA at 1:62500. ETV5 antibody can be used for detection of ETV5 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

### Expiration Date

12 months from date of receipt.



Antibody used in WB on Human 293T at 0.2...



Antibody used in IHC on Human Spleen at ...