

Product Datasheet

DLX3 Antibody (orb1243897)

Description

DLX3 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human DLX3.

Target

DLX3

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store DLX3 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

DLX3 antibody can be used for detection of DLX3 by ELISA at 1:1562500. DLX3 antibody can be used for detection of DLX3 by western blot at 1.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

32 kDa

Uniprot ID

[O60479](#)

NCBI


[NP_005211](#)

Dilution Range

DLX3 antibody can be used for detection of DLX3 by ELISA at 1:1562500. DLX3 antibody can be used for detection of DLX3 by western blot at 1.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date


12 months from date of receipt.



Antibody used in IHC on Human kidney at ...



Antibody used in WB on Transfected 293T ...



Antibody used in IHC on Human Placenta.