



Product Datasheet

TLR2 antibody (orb12439)


Description	Rabbit polyclonal antibody to TLR2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	A synthetic peptide of human TLR2(NP_003255.2).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	105KDa
Uniprot ID	O60603
Dilution Range	WB: 1:200-1:500, IHC-P: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of mouse lung tissue.



Immunofluorescence analysis of HeLa cells.



Immunohistochemical staining of mouse tissue.