

Product Datasheet

MYF6 Antibody (orb1243913)

Description

MYF6 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human MYF6.

Target

MYF6

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store MYF6 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

MYF6 antibody can be used for detection of MYF6 by ELISA at 1:1562500. MYF6 antibody can be used for detection of MYF6 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

27 kDa

Uniprot ID

[P23409](#)

NCBI


[NP_002460](#)

Dilution Range


MYF6 antibody can be used for detection of MYF6 by ELISA at 1:1562500. MYF6 antibody can be used for detection of MYF6 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date


12 months from date of receipt.



Antibody used in IHC on Human Pancreas.



Antibody used in IHC on Human Intestine.



Antibody used in WB on Human heart at 0....