

Product Datasheet

ZBTB48 Antibody (orb1243967)

Description

ZBTB48 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on r-ZBTB48
at 0.2-1...

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ZBTB48.

Target

ZBTB48

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store ZBTB48 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Antibody
used in WB
on Human
293T at 1 u...

Note

For research use only

Application notes

ZBTB48 antibody can be used for detection of ZBTB48 by ELISA at 1:1562500. ZBTB48 antibody can be used for detection of ZBTB48 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

77 kDa

Uniprot ID

[P10074](#)

NCBI

[NP_005332](#)

Dilution Range

ZBTB48 antibody can be used for detection of ZBTB48 by ELISA at 1:1562500. ZBTB48 antibody can be used for detection of ZBTB48 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

17 kDa

46 kDa

36 kDa

21 kDa

ZBTB48

80

65

45

31

22

Protein A/G-ZBTB48
Sample Type: Human 293T
Antibody Concentration: 1 µg/mL