

Product Datasheet

OR13C9 Antibody (orb1243979)

Description

OR13C9 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human OR13C9.

Target

OR13C9

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store OR13C9 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

OR13C9 antibody can be used for detection of OR13C9 by ELISA at 1:12500. OR13C9 antibody can be used for detection of OR13C9 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

36 kDa

Uniprot ID

[Q8NGS9](#)

NCBI


[NP_001001956](#)

Dilution Range


OR13C9 antibody can be used for detection of OR13C9 by ELISA at 1:12500. OR13C9 antibody can be used for detection of OR13C9 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in IHC on Human Brain at 4...



Antibody used in IHC on Human Kidney at ...



Antibody used in WB on Human HepG2 at 0....