

Product Datasheet

NKX2-8 Antibody (orb1243994)

Description

NKX2-8 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human NKX2-8.

Target

NKX2-8

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store NKX2-8 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

NKX2-8 antibody can be used for detection of NKX2-8 by ELISA at 1:12500. NKX2-8 antibody can be used for detection of NKX2-8 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

26 kDa

Uniprot ID

[O15522](#)

NCBI

[NP_055175](#)

Dilution Range

NKX2-8 antibody can be used for detection of NKX2-8 by ELISA at 1:12500. NKX2-8 antibody can be used for detection of NKX2-8 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

90kDa
80kDa

32kDa
25kDa

Antibody
used in WB
on Human
Heart lysat...

NKX2-8

90—
80—
65—
40—
31—
22—

Protein Ase-NKX2-8
Sample Type: Human Fetal Heart
Antibody Concentration: 1µg/mL

Antibody
used in WB
on Hum.
Fetal Heart ...

NKX2-8

90—
80—
65—
40—
31—
22—

Protein Ase-NKX2-8
Sample Type: Human Fetal Liver
Antibody Concentration: 1µg/mL

Antibody
used in WB
on Hum.
Fetal Liver ...