


Product Datasheet

GATA2 Antibody (orb1243995)


Description

GATA2 Antibody

Species/Host	Rabbit
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human GATA2.
Target	GATA2
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store GATA2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	GATA2 antibody can be used for detection of GATA2 by ELISA at 1:1562500. GATA2 antibody can be used for detection of GATA2 by western blot at 1.0-2.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	50 kDa
Uniprot ID	P23769
NCBI	NP_116027
Dilution Range	GATA2 antibody can be used for detection of GATA2 by ELISA at 1:1562500. GATA2 antibody can be used for detection of GATA2 by western blot at 1.0-2.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in IHC on Human Liver cell...



Antibody used in WB on Mouse liver, N2a ...



Antibody used in WB on Human K562 cells ...