

Product Datasheet

KIN Antibody (orb1244018)

Description

KIN Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human KIN.

Target

KIN

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store KIN antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

KIN antibody can be used for detection of KIN by ELISA at 1:2500. KIN antibody can be used for detection of KIN by western blot at 0.2-5.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

45 kDa

Uniprot ID

[O60870](#)

NCBI


[NP_036443](#)

Dilution Range


KIN antibody can be used for detection of KIN by ELISA at 1:2500. KIN antibody can be used for detection of KIN by western blot at 0.2-5.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date


12 months from date of receipt.



Antibody used in IHC on Human Spermatozoph...



Antibody used in IHC on Human kidney lys...



Antibody used in WB on Human HepG2 at 0....