

Product Datasheet

KLF8 Antibody (orb1244046)

Description

KLF8 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human KLF8.

Target

KLF8

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store KLF8 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

KLF8 antibody can be used for detection of KLF8 by ELISA at 1:2500. KLF8 antibody can be used for detection of KLF8 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

39 kDa

Uniprot ID

[O95600](#)

NCBI


[NP_009181](#)

Dilution Range


KLF8 antibody can be used for detection of KLF8 by ELISA at 1:2500. KLF8 antibody can be used for detection of KLF8 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in IHC on Human kidney lys...



Antibody used in WB on Human Kidney at 0...