

Product Datasheet

STON1-GTF2A1L Antibody (orb1244051)

Description

STON1-GTF2A1L Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Human
HepG2 at
5....

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SALF.

Target

STON1-GTF2A1L

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store SALF antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

SALF antibody can be used for detection of SALF by ELISA at 1:62500. SALF antibody can be used for detection of SALF by western blot at 5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

132 kDa, 52 kDa, 50 kDa

Uniprot ID

[Q9Y6Q2](#)

NCBI

[NP_758515](#)

Dilution Range

SALF antibody can be used for detection of SALF by ELISA at 1:62500. SALF antibody can be used for detection of SALF by western blot at 5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.