

Product Datasheet

YAF2 Antibody (orb1244061)

Description

YAF2 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human YAF2.

Target

YAF2

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store YAF2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

YAF2 antibody can be used for detection of YAF2 by ELISA at 1:62500. YAF2 antibody can be used for detection of YAF2 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

20 kDa

Uniprot ID

[Q8IY57](#)

NCBI

[NP_005739](#)

Dilution Range


YAF2 antibody can be used for detection of YAF2 by ELISA at 1:62500. YAF2 antibody can be used for detection of YAF2 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in WB on Human Jurkat 1.25...



Antibody used in IHC on Human Muscle.