

Product Datasheet

GTF2H3 Antibody (orb1244067)

Description

GTF2H3 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human GTF2H3.

Target

GTF2H3

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store GTF2H3 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

GTF2H3 antibody can be used for detection of GTF2H3 by ELISA at 1:2500. GTF2H3 antibody can be used for detection of GTF2H3 by western blot at 2.0/mL.

Clonality

Polyclonal

MW

34 kDa

Uniprot ID

[Q13889](#)

NCBI


[NP_001507](#)

Dilution Range

GTF2H3 antibody can be used for detection of GTF2H3 by ELISA at 1:2500. GTF2H3 antibody can be used for detection of GTF2H3 by western blot at 2.0/mL.

Expiration Date

12 months from date of receipt.



Antibody used in IHC on Human Muscle.



Antibody used in WB on Human Jurkat 0.2-...