


Product Datasheet

SPIB Antibody (orb1244074)

Description

SPIB Antibody

Species/Host	Rabbit
Reactivity	Canine, Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SPIB.
Target	SPIB
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store SPIB antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	SPIB antibody can be used for detection of SPIB by ELISA at 1:62500. SPIB antibody can be used for detection of SPIB by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	19 kDa
Uniprot ID	Q01892
NCBI	NP_003112
Dilution Range	SPIB antibody can be used for detection of SPIB by ELISA at 1:62500. SPIB antibody can be used for detection of SPIB by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in WB on Human Daudi at 0....