

Product Datasheet

HOXD4 Antibody (orb1244082)

Description

HOXD4 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on
Transfected
293T ...

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human HOXD4.

Target

HOXD4

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store HOXD4 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

HOXD4 antibody can be used for detection of HOXD4 by ELISA at 1:62500. HOXD4 antibody can be used for detection of HOXD4 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

28 kDa

Uniprot ID

[P09016](#)

NCBI

[NP_055436](#)

Dilution Range

HOXD4 antibody can be used for detection of HOXD4 by ELISA at 1:62500. HOXD4 antibody can be used for detection of HOXD4 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.