


Product Datasheet

KLF9 Antibody (orb1244087)


Description

KLF9 Antibody


Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human KLF9.
Target	KLF9
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store KLF9 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	KLF9 antibody can be used for detection of KLF9 by ELISA at 1:62500. KLF9 antibody can be used for detection of KLF9 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	27 kDa
Uniprot ID	Q13886
NCBI	NP_001197
Dilution Range	KLF9 antibody can be used for detection of KLF9 by ELISA at 1:62500. KLF9 antibody can be used for detection of KLF9 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in IHC on Human kidney at ...



Antibody used in WB on Human Fetal Liver...



Antibody used in IHC on Human Testis.