




Product Datasheet

POU3F2 Antibody (orb1244094)

Description	POU3F2 Antibody	POU3F2
Species/Host	Rabbit	 See IHC 2 Data and Customer Feedback for more information Antibody used in IHC on Mouse adult cort...
Reactivity	Canine, Human, Mouse, Rat	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human POU3F2.	
Target	POU3F2	 Antibody used in WB on Human Jurkat 5 ug...
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store POU3F2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	 Antibody used in IHC on Human Melenoma a...
Application notes	POU3F2 antibody can be used for detection of POU3F2 by ELISA at 1:12500. POU3F2 antibody can be used for detection of POU3F2 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	47 kDa	
Uniprot ID	P20265	
NCBI	NP_005595	
Dilution Range	POU3F2 antibody can be used for detection of POU3F2 by ELISA at 1:12500. POU3F2 antibody can be used for detection of POU3F2 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	