



Product Datasheet

GFI1B Antibody (orb1244124)

Description	GFI1B Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC, WB	Antibody used in IHC on Human kidney at ...
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human GFI1B.	
Target	GFI1B	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	Antibody used in IHC on Human Heart at 4...
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store GFI1B antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	GFI1B antibody can be used for detection of GFI1B by ELISA at 1:312500. GFI1B antibody can be used for detection of GFI1B by western blot at 1.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	38 kDa	Antibody used in WB on Human Jurkat 0.2-...
Uniprot ID	Q6T888	
NCBI	NP_004179	
Dilution Range	GFI1B antibody can be used for detection of GFI1B by ELISA at 1:312500. GFI1B antibody can be used for detection of GFI1B by western blot at 1.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	