

Product Datasheet

TFAP2D Antibody (orb1244126)

Description

TFAP2D Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat, Zebrafish

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Human
Muscle at 2...

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human TFAP2D.

Target

TFAP2D

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store TFAP2D antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

TFAP2D antibody can be used for detection of TFAP2D by ELISA at 1:312500. TFAP2D antibody can be used for detection of TFAP2D by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

50 kDa

Uniprot ID

[Q7Z6R9](#)

NCBI

[NP_758438](#)

Dilution Range

TFAP2D antibody can be used for detection of TFAP2D by ELISA at 1:312500. TFAP2D antibody can be used for detection of TFAP2D by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.