

## Product Datasheet

### ZBTB45 Antibody (orb1244141)

## Description

ZBTB45 Antibody

### Species/Host

Rabbit

### Reactivity

Human

### Conjugation

Unconjugated

### Tested

ELISA, WB

### Applications

Antibody  
used in WB  
on Human  
Lung at 0.0...

### Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ZNF499.

### Target

ZBTB45

### Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

### Form/Appearance

Liquid

### Concentration

batch dependent

### Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store ZNF499 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

### Note

For research use only

### Application notes

ZNF499 antibody can be used for detection of ZNF499 by ELISA at 1:312500. ZNF499 antibody can be used for detection of ZNF499 by western blot at 0.015 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

### Clonality

Polyclonal

### MW

54 kDa

### Uniprot ID

[Q96K62](#)

### NCBI

[NP\\_116181](#)

### Dilution Range

ZNF499 antibody can be used for detection of ZNF499 by ELISA at 1:312500. ZNF499 antibody can be used for detection of ZNF499 by western blot at 0.015 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

### Expiration Date

12 months from date of receipt.