

Product Datasheet

NFIC Antibody (orb1244152)

Description

NFIC Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human NFIC.

Target

NFIC

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store NFIC antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

NFIC antibody can be used for detection of NFIC by ELISA at 1:1562500. NFIC antibody can be used for detection of NFIC by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

55 kDa

Uniprot ID

[P08651](#)

NCBI


[NP_995315](#)

Dilution Range

NFIC antibody can be used for detection of NFIC by ELISA at 1:1562500. NFIC antibody can be used for detection of NFIC by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date


12 months from date of receipt.



Antibody used in IHC on Human Adipocytes...



Antibody used in WB on Human MCF-7 at 1 ...



Antibody used in IHC on Human Adipocytes...