

Product Datasheet

KIF20A Antibody (orb1244184)

Description

KIF20A Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human KIF20A.

Target

KIF20A

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store KIF20A antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

KIF20A antibody can be used for detection of KIF20A by ELISA at 1:1562500. KIF20A antibody can be used for detection of KIF20A by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

100 kDa

Uniprot ID

[O95235](#)

NCBI


[NP_005724](#)

Dilution Range


KIF20A antibody can be used for detection of KIF20A by ELISA at 1:1562500. KIF20A antibody can be used for detection of KIF20A by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date


12 months from date of receipt.



Antibody used in IHC on Human Thyroid at...



Antibody used in WB on Human 721_B at 1 ...



Antibody used in IHC on Human Uterus.