

Product Datasheet

LYPLA2 Antibody (orb1244206)

Description

LYPLA2 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human LYPLA2.

Antibody used in WB on Human 293T at 0.2...

Target

LYPLA2

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store LYPLA2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Antibody used in IHC on Human Adult Live...

Note

For research use only

Application notes

LYPLA2 antibody can be used for detection of LYPLA2 by ELISA at 1:62500. LYPLA2 antibody can be used for detection of LYPLA2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

25 kDa

Uniprot ID

[O95372](#)

NCBI

[NP_009191](#)

Dilution Range

LYPLA2 antibody can be used for detection of LYPLA2 by ELISA at 1:62500. LYPLA2 antibody can be used for detection of LYPLA2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

