


Product Datasheet

ERLIN2 Antibody (orb1244237)


Description

ERLIN2 Antibody


Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ERLIN2.
Target	ERLIN2
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store ERLIN2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	ERLIN2 antibody can be used for detection of ERLIN2 by ELISA at 1:62500. ERLIN2 antibody can be used for detection of ERLIN2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	38 kDa
Uniprot ID	O94905
NCBI	NP_009106
Dilution Range	ERLIN2 antibody can be used for detection of ERLIN2 by ELISA at 1:62500. ERLIN2 antibody can be used for detection of ERLIN2 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in WB on Human 293T at 0.2...



Antibody used in WB on Human HeLa, AT3 a...



Antibody used in WB on Hum. Adult Placen...