

Product Datasheet

CPE Antibody (orb1244240)

Description

CPE Antibody

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human CPE.
Target	CPE
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store CPE antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	CPE antibody can be used for detection of CPE by ELISA at 1:312500. CPE antibody can be used for detection of CPE by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	52 kDa
Uniprot ID	P16870
NCBI	NP_001864
Dilution Range	CPE antibody can be used for detection of CPE by ELISA at 1:312500. CPE antibody can be used for detection of CPE by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.


90 kDa

40 kDa

31 kDa

22 kDa

Antibody used in WB on Human MCF7 at 0.2...



Antibody used in IHC on Human Brain at 5...

CPE

90

65

40

31

22

Rabbit Anti CPE
Sample Type: Human Fetal Lung
Antibody Concentration: 1µg/mL

Antibody used in WB on Hum. Fetal Lung a...