

Product Datasheet

TRIM15 Antibody (orb1244269)

Description

TRIM15 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human TRIM15.

Target

TRIM15

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store TRIM15 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

TRIM15 antibody can be used for detection of TRIM15 by ELISA at 1:312500. TRIM15 antibody can be used for detection of TRIM15 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

52 kDa

Uniprot ID

[Q9C019](#)

NCBI

[NP_150232](#)


Dilution Range

TRIM15 antibody can be used for detection of TRIM15 by ELISA at 1:312500. TRIM15 antibody can be used for detection of TRIM15 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date


12 months from date of receipt.

TRIM15




Antibody used in WB on Hum. Adult Placen...

TRIM15



Antibody used in WB on Hum. Fetal Muscle...



Antibody used in WB on Human 721_B cells...